

- > 1 x Cardio Development Kit .
- >
- > Contents include: 3 x Main Shimmer Wireless Sensor Platforms
- >
- > 3 x ECG Modules
- >
- > 12 x ECG/EMG Leads
- >
- > 3 x Wearable Chest Straps
- >
- > 2 x USB Readers (This allows for programming and charging of the units)
- >
- > 3 x MicroSD cards (For storing captured data locally if desired)
- >
- > 1 x Manual
- >
- > 1 x Live USB (provides supporting applications, and quick start environment)
- >
- >
- >
- > 1 x Main Shimmer Unit .
- >
- > 1 x EMG Module .
- >
- > 1 x ECG/EMG Leads .

Simple configuration of three sensors of ECG, one sensor of Temperature and one sensor of SpO2

- >
- > With receiver of PC with your software.
- >
- > We do not have a SpO2 sensor. We are planning to release a GPS sensor by the end of January which will include a Temperature sensor in its functionality.
- >
- > To provide three Shimmer ECG units, the Cardio Development Kit should meet your requirements.
- >
- > The Development Kits are simply a combination of products, and are designed to meet a specific sensing requirements. For example, the Cardio Development Kit provides all the components necessary to set up an application that wishes to capture and transfer ECG readings wirelessly.

- >
- >